

Settings for your insulin (DIA and time-to-peak)

Contribution to the discussion among DIY loopers

The author assumes no liability

V.3.0 Sep'24



-
1. Setting insulin related parameters
 - 1.1 Mathematical model used
 - 1.2 Time-to-peak and DIA
 - 1.2.1 Insulin choice matters for profile ISF, IC
 - 1.2.2 Duration of insulin action
 - 1.2.3 Quantitative effects of changing DIA
 - 1.2.4 Shorter time-to-peak allows better control
 2. Other factors of potential relevance
 - 2.1 Age (of the diabetic)
 - 2.2 Dose
 - 2.3 Scatter (imprecision)
 3. Mixes of two insulins
 4. U200 insulin
 5. Utilization of insulins with super fast bio-availability
 - 5.1 i.v. insulin utilization
 - 5.2 Inhaled insulin (Afrezza)
 - 5.3 Lyumjev microGlucagon mix
 - 5.4 Development of novel super fast insulins
-

Before doing any other tuning, make sure you are on the insulin you really want to be on, and have reasonably set the insulin-related parameters for your looping system.

In case you are just starting to loop and need to “household” with your time, all you need from this paper should be just two messages:

- Select your insulin in AAPS configuration, and refer to the data given in [section 1.2.2](#) regarding DIA. To set it on 7 h is a fair guess for a start, if you are uncertain.
- You probably heard that Lyumjev or Fiasp are in principle *) best for looping.

*) from an activity kinetics standpoint (see [section 1.2.4](#))

In [section 1.2.1](#) it is explained why it is a good idea to switch at the start of your looping journey, rather than at some later point.

Changes between insulins with similar time-to-peak, like Fiasp -> Lyumjev, will be easier, and will not require much of an adjustment as in the example given in [section 1.2.1](#).

But of course you can switch at any later time, as well. Many prefer actually to start looping with a less reactive insulin. This enhances safety in the initial months of

getting to know, and tuning, the loop. (Same thought is behind the Objectives in AAPS that give access to SMBs only after a couple of other steps).

Also, struggling with too many occlusions (and pain) can make it difficult to switch to one of the fastest insulins.

1. Setting insulin related parameters

Besides time-(minutes) to-peak activity, also the duration of insulin action (DIA, hours) that you select in your profile strongly influences how the loop calculates the activity from insulin, as it unfolds in every 5-minute segment that your loop analyzes.

1.1 Mathematical model used

Especially what should be selected as duration of insulin action (DIA) is very strongly influenced by the model used to figure out active insulin two, three, and more hours after administration. Misunderstandings about this is often a source for disputes between loopers and their treating physicians.

All insulin administrations (bigger and minor) add up to a insulin activity pattern. In the case of looping, with user boli, basal insulin, TBR modifications and SMBs given at various times, with overlapping DIAs, this can be quite complex.

In AAPS you can see insulin activity in your main screen as an extra thin yellow curve. Together with carb absorption is “explains” most of what you see in your glucose curve.

This **insulin activity pattern is an extremely important basis for each of your loop’s decisions**. Having the wrong settings would give your semi-automated insulin management a permanent drift towards over- or towards under-corrections.

The loop system can still counter-regulate, but – if you burden your’s with wrong DIA or time-to-peak settings in your profile – this would “use up” some of it’s (limited) capacity to regulate for you.

Example: After heavy dinner, a DIA set too short “tells your loop” that active insulin is practically gone after time X. The loop takes that info for granted, and if it sees some insulin needed at that time X (and be it only for your profile basal need - as you also

communicate to the loop, you need to remain stable -), then, at night-time, the loop will give you more insulin than you really need.

Therefore, before you tune your ISF differently, make sure to have a look at your DIA setting.

Please understand (and see to it, that your treating professionals understand) that models can differ strongly:

- DIY looping systems use the – less common – exponential decay model.
- Medtronic uses non-linear capped curves (as in handbook to their pumps)
- Doctors / diabetes educators mostly have a rough linear model in mind
- xDrip uses a bilinear math (“with kinks”) to model insulin activity (Caution: This info might be outdated)

All models are working „good enough“ for their (main) intended applications. But, as explained above, it is worth the effort to use an exact modelling of insulin activity for a loop, so it can perform optimally.

The mathematical model of insulin activity over time that we use anchors on time-to-peak (minutes) and on DIA (hours) in characteristic ways. This is quantitatively shown for exponential decay models in [section 1.2.3-](#)

In AAPS, the insulin tab shows two curves:

The **pink** one starts at 1.0 (100%) and goes down to 0 (0%) when the DIA is over. It shows how much of the total activity (the capacity to lower bg) is left, at any time. So, it is like the job number we always have in our AAPS home screen. The problem with that, as with the pink curve, is that it may give you a false impression regarding how much “power” there actually is, now, as you need it. That is where the other curve (and on your AAPS home screen, the related thin yellow insulin activity curve) come in:

The **blue** one shows how the activity goes: Practically nothing (!) for a bunch of minutes, then rapidly going high, and then slowly fading out over the DIA period (with a maximum at time-to-peak). For its calculations, AAPS adds these blue curves up for all boli, SMBs and TBRs ≠ profile basal!


1.2 Time-to-peak activity and DIA for various insulins

Principally, there are „correct“ settings specific for each insulin type, notably regarding time-to-peak activity. This is pre-programmed in the insulin choices for AAPS, for instance.

Regarding the DIA to set, there is more uncertainty. The following mostly cites or summarizes findings published by David Burren.

1.2.1 Insulin choice matters for profile ISF, IC

The following chart is *the inverse* of the pink curve in the AAPS insulin tab: *Not insulin still there to be used*, but Insulin used up, going from 0% towards 100% in the 9 h DIA, for Humalog with 75 minutes, and for Lyumjev with 45 minutes time-to-peak.



From a simplistic point of view, you can see that at the two-hour mark, more of the Lyumjev (75.5%) should have had effect than the Humalog (50.2%).

So when we're calculating *how much insulin to give for a correction*, we should tell it to give more Humalog up-front to get the same result after 2 hours.

“With the current incarnation of the oref1 algorithm I do find it helpful to scale the ISF (and IC) relative to the amount of each insulin used in the first 2 hours. For Fiasp ~66% is used in the first 2 hours. Lyumjev 75% (and NovoRapid/etc 50%).

Thus for Lyumjev I use ratios that are 15% larger than for Fiasp.

Basal rates are unchanged.” (D. Burren, AAPS Users 03Apr.2024)

The system will of course be tracking the IOB and forecasting the BG curves for hours into the future, so we do have some safety built in regarding the extra insulin.

When going from Humalog to using Lyumjev, this must have some consequences for the Insulin Sensitivity Factor (ISF) to use in the profile. If, for example, you had 1.8 mmol(I/U for

123 Humalog, you should expect a “good ISF for going with Lyumjev” in the area of 2.7 mmol/l/U.
124 According to the curves shown above (at dotted 2 hr line) a factor 75.5/50.2 applied yields
125 the same amount of insulin for a correction.

126 Likewise, the Carb Ratio (IC) may deserve an adjustment when switching insulins.

- 127 • The IC could be adjusted by the same factor, for instance it might go from 7.7 g/U
128 (Humalog) to 11.6 g/U (Lyumjev).

129 For a meal of 60 g, 7.8 U (=60/7.7 g/U) Humalog would have contributed 3.9 U
130 (=50.2%*7.8 U); likewise, 5.2 U (=60g/11.6 g/U) would have contributed 3.9 U
131 (=75.5%*5.2 U)

- 132 • For meals bigger than about 60 g you should observe that, while your insulin bolus has good
133 activity, only a limited number of carbs can get digested (30 g/h seems the limit for most).
134 Refer to the paper on IC determination (section “Determination at meal times”) in:
135 <https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-settings/tree/FCL-w/autoISF>

136 The given example showed that switching to a “faster” insulin can have relevant
137 consequences for your key profile parameters.

138 David Burren also reports that between the two rather extreme insulin choices he tested, the
139 total amounts of insulin (TDD) did not significantly differ (- as we would expect: The same
140 amounts just gets delivered slower, even at same selected DIA, with Humalog).

141 But while the TDD has *not* changed, the instantaneous insulin levels *have*.

142 When the system is fighting post-meal “highs” the IOB will be noticeably lower with Lyumjev.

143 Although the average overall level remains similar, this might have some implications
144 for the concept of hyper-insulinaemia. This may be a subtle advantage of faster
145 insulins.

146

147 1.2.2 Duration of insulin action

148

149 The following focusses on the more uncertain topic of which duration of insulin action (DIA)
150 to use. It largely relies on, and quotes, results from several thorough investigations done
151 by David Burren: (<https://bionicwookiee.com>)

152


153 The numbers he ended up preferring are:

Insulin	Peak	Duration (DIA)
Humalog	75 minutes "Rapid-Acting Oref"	7 hours *)=
NovoRapid	75 minutes "Rapid-Acting Oref"	9 hours
Fiasp	55 minutes "Ultra-Rapid Oref"	9 hours

154 *) Later investigations <https://bionicwookiee.com/2022/04/13/revised-humalog-model-in-a-closed-loop/>
155 led to suggesting 9 hours DIA also for Humalog

156 The default constraints in AAPS have the duration limited to 7 hours, so he had to
157 make some local changes to the limits. It's also possible if you set your "patient type"
158 to "Pregnant", but if so you need to carefully check all the affected safety limits ([listed](#)
159 [in the AAPS documentation](#)). This may change in a future update to AndroidAPS.

160 For Lyumjev (45 minutes; Lyumjec Oref), there is not a big difference between a 7 and a 9 h
161 DIA:



162 However, David Burren ([https://bionicwookiee.com/2022/04/13/revised-humalog-](https://bionicwookiee.com/2022/04/13/revised-humalog-model-in-a-closed-loop/)
163 [model-in-a-closed-loop/](#)) reports that, despite it's a very subtle change, he has found
164 it can make a significant difference 5-6 hours after a meal, when ...the tails of the
165 earlier doses do add up, and the system had been underestimating the IOB when
166

calculating (using the shorter DIA) what was needed with new doses. With changing to a longer DIA, his average Time Below Range has reduced. (Comment: This is interesting for “fine-tuning extremists” but probably only a formal gain of little clinical significance, assuming the bg curve just swinging a bit more often, or longer, by a few mg/dl below the 70 mg/dl cut-off, that defines “below range”. Judge from your own data, when/if fine tuning.)

On the DIA topic for various insulins see also: <https://www.diabettech.com/insulin/why-we-are-regularly-wrong-in-the-duration-of-insulin-action-dia-times-we-use-and-why-it-matters/>

1.2.3 Quantitative effects of changing DIA

Any given insulin dose comes with a defined total capacity for a certain bg lowering effect. How strong or weak this unfolds over a couple of hours can be mathematically modelled. Inoref(1) systems, time-to-peak and DIA completely define this curve.


We can look on effects of increasing the set DIA in terms of how insulin activity would differ at any moment after administering a dose.

The next example given (chart below) does that for going from a 5 h DIA, a 6 h DIA or a 7 h DIA towards 9 h for Lyumjev

We see the peak going lower, and the tail activity higher when DIA is increased:

LYUMJEV	peak @45m	max effect on “tail” at ~ 3.5 h (220 minutes)
DIA 5 → 9h	minus 5.5 %	plus 2.5%
DIA 5 → 6h	minus 2,7 %	plus 1.3%
DIA 6 → 9h	minus 2.8 %	plus 1.2%
DIA 6 → 7h	minus 1,4 %	plus 0.6%
DIA 7 → 9h -	minus 1.4 %	plus 0.6%

So, the „tail“ effects differ by less than 3 percent (of peak activity=100%) in the later stages of DIA:



While 3 % sounds low, the significance of the problem should not be underestimated:

- For our Lyumjev case, note that the quoted 3% result is 3% of maximal activity.
Example: Activity at 180 minutes is about 0.0010 compared to 0.0080 at peak (blue curve in AAPS INS tab). 2.5% of 0.0080 would be 0.0002. BUT: 0.0012 is 20 % more than 0.0010, so **REALLY the difference in activity at 180 minutes can be up to 20%** . Still, after a bolus of 8 units (and/or SMBs that reach that iob level) for a typical meal, the max. difference from 5 -> 9 hour DIA would roughly be, whether 1.0 U or 1.2 U are active iob left at 180 minutes. That **difference (+ 0.2 U)** should be within the loop's regulating capacity from reducing basal.
 However, it becomes much bigger for users of other insulins (with longer time-to-peak):
- The delta **effects get much bigger with insulins that have a longer time-to-peak**
 Some quantitative data for other insulins are as follows:

FIASP (peak=60m) min/max differences

DIA 5 → 9h | 6 → 9h | 7 → 9h: -10,1 / **+6,8%** | -5,6 / +3,0% | -2,9 / +1,4%

NOVORAPID (peak=75m) min/max differences

DIA 5 → 9h | 6 → 9h | 7 → 9h: -15,4 / **+14,1%** | -9,1 / +7,0% | -4,8 / +3,0%

Above example applied to Novorapid **): The effect would be up to +14.1% of max (!)
 => 2.1 U instead 1 U at 180 minutes. A **difference of + 1.2 U** results here, if DIA is
 set at 5, not at 9 h, so **REALLY** it could go **up to + 120% !**)

More see: szantos, de.loopercommunity.org May 2022

<https://de.loopercommunity.org/t/naechtlicher-unterzucker/10626>

**) 2,5% → + 0.2 U ergo 14.1% → +1,1 U stimmt insofern nicht ganz genau, als man beim Novorapid Case auch die Novorapid Peak-Höhe zugrunde legen müsste (die ich aber nicht greifbar habe). Wenn diese von Haus aus 20% niedriger nur kommt, hätten wir ca + 0.9U, also weiterhin etwa eine Verdoppelung ... die wir mit unserer Wahl eines längerem DIA unserem Loop sagen könnten, damit er entsprechend weniger zu-schiesst ... ergo weniger Hypogefahr hinten heraus ...

1.2.4 Shorter time-to-peak allows better control of meal-related bg spikes

A modelling study „The Artificial Pancreas and Meal Control“ by A. El Fathi et al

([https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-settings/blob/HCL-.settings-main-repo-\(pdf\)/The%20Artificial%20Pancreas%20and%20Meal%20Control.pdf](https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-settings/blob/HCL-.settings-main-repo-(pdf)/The%20Artificial%20Pancreas%20and%20Meal%20Control.pdf)) can help us understand the effects on glucose peak heights from the course of carb absorption and of insulin activity. The graph shows on the y axis peak over baseline (the overall deltaBG in mmol/l), and on the x-axis the relative speed of insulin absorption to carb absorption. Carb absorption is always faster, therefore all values are under 1.0. But with Lyumjev we move closer to 1. The model calculation shows that **faster**

insulins (red dotted) will result in **lower** glucose **peaks** than slower insulins (violet dotted: reduction by 46% or minus 2.5 mmol/l = minus 46 mg/dl after a 60 g carb load).

So, this model supports that using a faster insulin will

- lead to less high glucose peaks, notably for bigger meals
- or might **tolerate** a couple of minutes **delayed** first meal bolus while not incurring unacceptable height of peaks.




FIGURE S2 A plot of the maximum glucose peak after ingestion of different carbohydrate (CHO) quantities as a function of the ratio between time-to-peak of meal absorption τ_m and time-to-peak of insulin absorption τ_i . This graph shows that, for instance, following a 60-g meal, the maximum peak of glucose is 5.4 mmol/L for a ratio $\alpha = \tau_m/\tau_i = 0.8$. Increasing the ratio to 0.9 (by slowing the meal digestion or providing a faster-acting insulin) may result in decreasing the peak by 46% to 2.5 mmol/L.

The latter is a pre-requisite for full closed loop, in which we leave it up to the loop to notice that a meal „must have started“, and to come forward with SMBs that are typically delayed compared to the bolus as given in hybrid closed loop. This is an encouraging result.

Moreover, the same chart shows us that the spread between the colored curves (they stand for different meal sizes) becomes significantly smaller when we move to the faster insulin with a 0.9 ratio. This means the danger of increasingly high post-meal glucose peaks for high-carb meals is sharply reduced, too. For example, the green curve suggests with the „0.8 insuline“ a peak of 10mmol/l (180 mg/dl) above your glucose level at meal start, but only +4 mmol/l (+70 mg/dl) when using a faster insulin with factor 0.9, which, when starting at or under 110 mg/dl, could keep glucose in range.

The message we can take from this is:: **The higher carb loaded our diet, the more important to use the fastest-available insulin.**

Note that this section 1.2.4 was about *“the real”* time-to-peak coming with your insulium choice; while the preceeding sections on DIA were about effects if *your set* DIA is more or less “off” reality.

2. Other factors of potential relevance


Source: [szantos](#)

The findings reported below can give you hints in which direction to look if you attempt to fine-tune your settings further, from the standard suggestion what should be suitable for your insulin ([section 1.2.2.](#)).

2.1 Age (of the diabetic)


ema.europa.eu

[novorapid-h-c-258-p46-0044-epar-assessment-report_en.pdf 3](#)



2.2 Dose

<https://journals.sagepub.com/doi/10.1177/1932296813514319>



[image1000×632 78.4 KB](#)

2.3 Scatter (imprecision)

Individual deviations from standard suggestions could also be justified by the fact, that all studies that underly the previously reported suggestions, come with very significant person-to-person scatter.

All lines in the charts, as above shown from studies, are averaged data. (Some studies are indicating the very significant scatter seen, as well).

https://www.researchgate.net/figure/Pharmacodynamic-profiles-Insulin-action-as-expressed-as-GIR-required-to-maintain_fig1_41424712

3. Mixes of two insulins

The author did for some time successfully use a 50/50 mix of Fiasp and Novorapid, applying the time-to-peak for Fiasp, and longest of the two DIA, as was suggested at the time, for these insulins.

For a more thorough discussion see <https://bionicwookiee.com/2022/03/02/mixing-insulins-theory-and-practice/>

and also: <https://bionicwookiee.com/2023/06/03/arcane-lyumjev-experiments/>

4. U200 insulins

Using up-concentrated insulins, e.g. in a U200 form, is sometimes chosen by loopers

- to reduce needed daily insulin volumes and get more time from 1 pump filling (pod)
- to reduce volume per injection for getting better tolerance regarding occlusions or pain

There are no relevant effects on insulin parameters like DIA and time-to-peak.

However, dilution or up-concentration factors are highly relevant for setting profile factors like ISF and IC, and also for some important safety settings like max iob for instance.

Refer to special discussions on that topic, e.g. here re. U200 Lyumjev
<https://www.diabettech.com/lyumjev/living-with-lyumjev-almost-a-year-in-review/> :

and also: <https://bionicwookiee.com/2023/06/03/arcane-lyumjev-experiments/>

5. Complementary utilization of insulins with super fast bio-availability

The effect of time to peak activity on bg control was shown quantitatively in the study presented in [section 2.1.4](#). It is the core problem coming with any sub-cutaneous insulin provision (via sub-cutaneous injections, or via insulin pumps), that time-to-peak activity can be rather long.

In the following we touch on faster ways to get insulin into blood. **Please note that the author does not encourage you to try any of those** unless in a medically supervised study context. Also, be aware that faster acting insulins further increase danger of hypoglycemia.

5.1 Using i.v. insulin

The “beauty” of using intra-venous insulin at around meal starts would stem from its much faster bio-availability, and also much shorter DIA

Case report: I do a basically unknown amount of insulin intravenously, let it be between 4 and 8 units (well below the size that my meal bolus would be). It really doesn't matter anyhow. What it does ,it brings me down to target within 30-40min. I record something like 4-6 units (so my loop doesn't want to get excessive insulin). Essentially, this prevents getting insulin longer than it actually has an effect (mine is gone from the system after 35min). To eliminate the false “activity tail” assigned also to the i.v. potion of insulin on bord, you can delete the i.v. insulin amount from the system *after it has done it's job* (not good for statistics/history data, but right, going forward without the DIA tail = letting your loop know the real iob).

It's an edge use (experimental) case . (source: Robert, discord FCL/iaAPS w autoISF, March 2024):

Please observe that this is not a recommendation to experiment with i.v. insulin unless in a medically supervised research context.

i.v. insulin is usually restriced to the surgical and intensive care hospital environments!

5.2 Inhaled insulin (Afrezza)

Afrezza is an inhaleable very fast (and also short) acting insulin which some find useful to correct high glucose levels.

Pro: An insulin inhalation can prevent, or quickly correct, a high bg, while avoiding a late hypo tendency from-a long “tail of effects” (DIA).

Cons: 1) Afrezza spray is hard to dose, and spraying into the lungs may not tolerated well. High cost. 2) It is not advisable to enter data into loop because the kinetics of this insulin are very different. => The short term problem is solved, but there will be consequences from skewed calculations if your looping mode relies on data like TDD, cob, Autosens, Autotune.

Conclusion

- 1) "Afrezza is fast enough to take effect *and also to get out (super short DIA)* that it doesn't necessarily impact the rest of the model: The spikes you'd normally expect don't tend to happen" (Tim Street).
You could look at it as, via Afrezza, "faking" "low carb meals always" to your loop: Your Afrezza shot takes quick care of the fast carbs (which you better not declare). All AAPS has to do, is take care of the effects from *slow carbs and FPU*s that need insulin when Afrezza effects are over.

(More discussion on Afrezza in AAPS looping see [Craig Gordon](#) 14Sep.2024 in Facebook AAPS users· [ernoSodstp6S40i 5heba14rf5m3hiip q i:1tm e708M6f7la0f2cl2Ptehi4](#)).

- 2) Primary approach should be to go without Afrezza, and avoid high bg by finding a proper meal management strategy (pre-bolus time, EatingSoonTT).
- 3) If you have access to Afrezza, try it occasionally *within your conventional meal management (2)* to see how it might help you manage rapid carbs.


Solving the problem at hand as best as we can, even if it makes the time afterwards a bit more complicated, is the name of the game. We and our loop do this all the time, for instance by giving more upfront insulin, then reducing basal (zero-temping). (quoted from slide 38 of: Meal Mgt. Basics, <https://github.com/bernie4375/HCL-Meal-Mgt.-ISF-and-IC-settings>)

Caution: Your looping algorithm is built on having realistic iob data, and carb (or carb deviation) data. **With Afrezza, you strongly interfere. It is your responsibility as early tester to think things through, so you can avoid risks.** For example:

- You probably must reduce the maxIOB set in AAPS Preferences by your usual (not-declared) Afrezza dose
- Timing of Afrezza application will make a huge difference in post-prandial bg curve characteristics. Applying as early as possible is super crucial if you operate with aggressive SMB settings (notably in autoISF FCL, where you probably must significantly lower the iobTH).
- But applying very early comes with immediate hypo risks. Be aware that, the faster acting the insulin, the more difficult to counter-act with glucose tablets!

5.3 Lyumjev + microGlucagon mix

A highly experimental approach to further improve the kinetics of Lyumjev for looping has been researched by a T1D physician (Dr. Carlsen, Trondheim, Norway) in self-experiments, with some promising first results.




Full article from Norwegian hospital researchers group:

<https://norwegianscitechnews.com/2024/01/a-new-insulin-mixture-could-make-life-easier-for-patients-with-type-1-diabetes/>

Our DIY researches pushed the topic further

(<https://discord.com/channels/953929437894803478/1253713103900639262>)

and the red curve may resemble their current status (graphs are scaled to show the relative strengths or patterns):



TrG+: Trondheim stronger G mix, tested on pigs

TrG: Trondheim weaker (?) G mix (equivalent to just lyumjev elsewhere.. their lyumjev control was way worse!)

Ly4: the 4th lyumjev article I had.. best approximation I think

GL: Trondheim original graph (dark blue here; was orange on previous page) which was the human one I think

GL50+M (red): modified GL graph of my own **DIY mix***) and blood tests for sequential insulin assays.

*) high dasiglucagon (200ng/unit) similar Trondheim high G , but other additives as described in Discord

Caution: Changing medications can be dangerous, and is not advised outside of a supervised medical study.

Also, your looping algorithm is built on having realistic iob data. By mixing-in glucagon, you strongly interfere. It is your responsibility as early tester to think things through, so you can avoid risks. For example, you should give AAPS a suitable insulin kinetic model (and make sure your mix is stable, regarding that, during the period of use).

See more related discussion:

- in Discord: <https://discord.gg/eHSgx5jWuk>
- in "AAPS Users" FB (Alf Einar Johnsen, Jan 19, 2024) <https://www.facebook.com/groups/AndroidAPSUsers/permalink/3733005573587499/?mibextid=W9rl1R>
- or in "iAPS unchained" FB (Robban Jansson, 06 May, 2024) <https://www.facebook.com/groups/151989761311250/permalink/286074914569400/>

5.4 Development of novel super-fast insulins

According to <https://arecor.com/wp-content/uploads/2022/01/220120-Areco-Formulation-Expertise-Exemplified-in-Diabetes.pdf> there is already a super fast insulin for sub-cutaneous pump delivery in clinical trial stage, AT247 (by Atera). It showed a 17 minutes **onset** of glucose-lowering action (compare: 37min.Novorapid, 23 min.Fiasp, data not given for Lyumjev, ca 20 min.). Unfortunately, it still has a long DIA and therefore it will not be the

424 ultimate thing (... and might turn out inferior to glucagon addition to Lyumjev as discussed in
425 [section 5.3](#))

426 So, that one could be a gradual improvement over Lyumjev.

427 Getting insulin activity from (micro)boluses something like 5 minutes earlier would further
428 limit bg rises (see [section 1.2.4](#)). However, it could not contribute much to resolving a major
429 problem we still see today in Full Closed Loop, which is: Fighting the rising bg is principally
430 limited because the inevitable “tail” of insulin activity produces hypo danger.

431 “Dreamspace”

432 In the following, as brief look follows into the **potential with a real break-through**, an even a
433 bit faster onset PLUS also a very much shortened tail of action.

434 An insulin that **acts even faster than any carb absorption** will, by design, be geared to
435 reacting on development of the bg curve. This brings huge **advantages**:

- 436
- No need for carb counting, and forbidden (!) to give meal boli.

437 Therefore such insulin should be restricted for pump use only.

- 438
- Rapid onset, PLUS short “tail” of activity, allows to operate near-physiologic data
439 regarding both, very small very frequent insulin doses, and bg not peaking super
440 high, ever.

441

- 442
- Technically this concept should be better than the dual hormone route (section 13.6 of
443 FCL e-book: [https://github.com/bernie4375/FCL-potential-autoISF-research/blob/FCL-book-](https://github.com/bernie4375/FCL-potential-autoISF-research/blob/FCL-book-autoISF)
444 [autoISF](https://github.com/bernie4375/FCL-potential-autoISF-research/blob/FCL-book-autoISF)) because it needs no additional reservoir and pump. Also, one reactive, fast
445 regulation is better than a sluggish regulation with counter-regulation (in form of a
446 quick correction by small titered dose of glucagon) on top.

447 However, the dual hormone system might be safer in real life – see first point under
448 ...problems (below).

449 But it also comes with a couple of **problems**:

- 450
- It requires to educate people that using that insulin in a syringe (or also manually
451 issuing a bolus from their pump) is an absolute no-go.
452 However, today everybody is used to doing just that, and we all know too well the
453 temptations to occasionally give a correction bolus.

454 It will be interesting to see how manufacturers would cover themselves from lawsuits
455 because normal people in some everyday situation, or in a special situation, may do it
456 in a way that now, with that innovative insulin, is definitely wrong. Note that once that
457 new insulin is “in”, no glucose tablet can rescue you any longer. You probably would
458 need glucagon syringe immediately.

459 (As already hinted on above, this is where the competing dual hormone approach to
460 FCL can play out an important advantage: It has that glucagon safeguard “right on
461 board”!)

- 462 • It absolutely requires a 24/7 reliable CGM (because the rapid corrections shoot
463 towards a target value that might be dangerously skewed)

464 Again, from everyday observations we know that only a minority even of today’s
465 loopers do have reliable CGM values every day. If a broader group of T1Ds is
466 targeted, this can turn out a real barrier.

- 467 • It also requires nearly 24/7 Bluetooth connectivity (because system will fall back to
468 profile basal, or to no insulin at all, when it cannot loop, and high bg would result,
469 coming with the temptation to issue a manual bolus....)

470 A pretty much hands-off (no bolus, no carb inputs) Full Closed Loop system could be
471 designed with a high degree of system integration and miniaturization
472 (pod+CGM+algo); this could provide rock-solid connectivity.

473 From a **commercial standpoint** this must make us wonder whether it would be commercially rolled
474 out at all, and for which target group.

475 It would be nice if it were launched to transform current (and future) commercial loops into FCL mode.
476 Due to liability concerns, especially in the US, very strong safety measures would be built-in, like not
477 allowing to set “low” target bg (limitations we already know from current commercial systems).
478 Commercially seen, this might open a big enough market, to offer the majority out there on high
479 HbA1c a rather care-free option to get into the ~ 7% range. If cost-effectiveness (and safety) can be
480 shown (and it might), this could revolutionize T1D treatment, with Full Closed Loop on a pump+CGM
481 becoming the standard. (The author is not familiar enough with type 2 diabetes management to say
482 anything whether, after some foolproof system integration (pod+CGM+algo), it could extend into that
483 very broad market, as well).

484 Especially if safety issues remain, alternatively routes to achieve around 80% TIR for a broad
485 T1D audience have already been described in a clinical study based on AAPS:
486 <https://pubmed.ncbi.nlm.nih.gov/36826996/>

487 Current **FCL users on autoISF** who like to see HbA1c around or even under 6% are a “too small
488 group” to be of any interest; but they could benefit from a super fast insulin (fast on-set PLUS fast
489 fading-out!) immediately, after giving their system the applicable kinetic curve for that insulin (in AAPS
490 easy possible under Config.Builder/Insulin/Free Peak Oref). This would improve fighting post-prandial
491 peak heights without ensuing hypo danger (which, currently with Lyumjev, is a delicate balancing act
492 in initial system tuning).